

OUTLINE

CAN ACCOMMODATE UNMOUNTED OPTICS:
 MINIMUM DIAMETER = Ø34.5mm
 MAXIMUM DIAMETER = Ø38mm
 MAXIMUM THICKNESS = 3.5mm

UNLESS OTHERWISE SPECIFIED
 DIMENSIONS ARE IN INCHES [MM]
 AND INCLUDE CHEMICALLY
 APPLIED OR PLATED FINISHES

SCALE: 2:1

WEIGHT: 0.06 lbmass

DR BY: J. WILLIAMSON


CH BY: DEROSA

PROJECT #: E3548

DATE DR: 12/14/2018

PROPRIETARY AND CONFIDENTIAL

THE INFORMATION CONTAINED IN THIS
 DRAWING IS THE SOLE PROPERTY OF
 NAVITAR. ANY REPRODUCTION IN PART
 OR AS A WHOLE WITHOUT THE WRITTEN
 PERMISSION OF NAVITAR IS PROHIBITED.



DO NOT SCALE DRAWING



RESOLV, FILTER HOLDER INTERFACE

SIZE	PART NUMBER	REV
A	1-81310	A